

Abstract

Capture and isolation of circulating metastatic cancer cell clusters (MCCC) is the essential first step to identify novel therapeutic targets focused on this important route of metastasis. We developed a unique biomimetic capture approach that allows selective isolation of these deadly clusters from the whole blood of lung cancer patients. Utilizing a microfluidic platform we developed a versatile smart-coating that combines biomimicry and immuno-capture for a highly efficient and selective isolation of viable MCCCs. Here we show the capture and characterization of multiple MCCCs from non-small cell lung cancer (NSCLC) patients. Table 1 shows the count of clusters detected in NSCLC patients based on a stringent criterion confirming their identification. The MCCCs were detected in every patient tested but none were detected in normal blood controls from donors of similar age. Isolated MCCC where characterized by immunofluorescence and submitted for genomic analyses.

Metastasis remains one of the leading causes of cancer related deaths worldwide. Insights into how the tumor propagates to pre-metastatic sites is necessary to develop focused therapeutic approaches. Circulating metastatic cancer cell clusters (MCCCs) have been established as a primary entity causing distal metastases. Significant advances have been made linking cluster integrity and viability to increased metastatic potential highlighting the unmet need for cluster focused therapies.

Our novel approach allows the rapid and easy capture of MCCCs paving the way for the discovery of novel or unique targets that eventually will bring new anti-metastatic therapies to the clinic.

What are MCCCs?




Figure 1: Basic concept of MCCC-driven metastasis. The primary tumor releases MCCCs and single circulating tumor cells, which travel through the bloodstream and localize elsewhere in the body forming potential metastatic seeds.

- MCCCs are also referred to as circulating tumor micro-emboli, CTC clusters, and collective cell migration.
- MCCCs migrate from the primary tumor as homotypic or heterotypic groups of cells intravasating into circulation.
- Contrasting to single cells and homotypic MCCCs, heterotypic MCCCs have been observed to preferentially escape immune-surveillance and survive in circulation, demonstrating a greater metastatic potential.
- MCCCs can get 'stuck' in micro-vasculature or actively extravasate into tissue forming pre-metastatic niches that may evolve into metastatic lesions.

Clinical significance of MCCCs




Figure 2: A brief history of circulating CTCs/MCCCs.

A snapshot of published clinical data correlating MCCCs with OS and PFS




Figure 3: Selected OS/PFS graphs from recent key publications on CTC-clusters/MCCCs detected in clinical patient studies.

The presence of MCCCs has a dramatic impact on the Overall Survival (OS) or Progression Free Survival (PFS) for cancer patients across most solid tumor types.

Sample processing





Figure 5: Schematic of sample processing. Syringe containing unprocessed patient blood sample is injected into MCCC capture chip at a flow rate of 100µL/minute. Effluent can be collected and stored for subsequent multiplexed analyses (ctDNA, exosome, etc.).



Capture optimization

- Co-cultured spheroids were generated by seeding A549 and HCC827 cells respectively, in combination with HFF-1 fibroblasts into a 384 well U-bottom plate.
- Spheroids were fluorescently stained and spiked into an 4mL freshly collected normal donor blood sample.
- The spiked blood sample was processed on a SMART-coated™ capture chip and washed with buffer.
- An accurate spiked spheroid count was taken and compared to the spheroids captured on the chip in order to calculate a capture efficiency percentage.
- Multiple normal blood donor samples without spiked clusters were processed and no clusters were observed.

Capture Efficiency using 1 mAb vs 3 mAb




Figure 7: Capture efficiency comparison using varying numbers of tissue specific capture antibodies. The unique properties of the dual-capture approach allows targeting of multiple antigens concurrently.

Capturing MCCCs from patients.

Whole blood samples of patients¹⁰ with confirmed lung cancers were collected and processed on a SMART-coated™ MCCC capture chip. Fast processing times keep MCCCs viable and naïve. Captured MCCCs were detected using automated bright field microscopy and either stained with immunofluorescence marker confirming MCCC-status or processed for genomic sequencing analysis.

Patient	Age/Gender	Diagnosis	Stage	# clusters detected	Sample use
INMC-044.1	82 / M	lung adenocarcinoma	IVb	2	MCCC detection by bright field microscopy and confirmation by immunofluorescence staining for CD44, EGFR, Hoechst
INMC-044.2	53 / F	squamous cell lung carcinoma	IV	1	
INMC-044.3	59 / F	metastatic lung adenocarcinoma	IV	2	
INMC-044.4	65 / M	squamous cell lung carcinoma	IV	1	
INMC-044.5	67 / M	non-small cell lung carcinoma	IIla	3	MCCC detection by bright field microscopy, samples processed for genomic analysis by RNA and DNA seq.
INMC-044.6	64 / F	lung adenocarcinoma	IV	MCCC were not IF stained to preserve RNA quality	
INMC-044.7	46 / F	non-small cell lung cancer	IV		
INMC-044.8	66 / F	non-small cell lung cancer	IV		

Table 1: Details of patient samples processed.

Immunofluorescence staining was performed *in situ* on the capture chip. This minimizes the sample handling and accelerates analysis. Imaging was performed on an automated Yokogawa CQ1high-content analysis system. Detected cell clusters were classified into three groups, [Hoechst positive, EGFR positive, and CD44 negative], [Hoechst positive, CD44 positive, and EGFR negative], and [Hoechst positive, CD44 positive, and EGFR positive]. Cells positive for only Hoechst and CD44 were presumed to be of hematopoietic origin and disregarded. CD44/EGFR staining without positive nuclear Hoechst stain were disregarded. Only clusters positive for all three stains were considered potential MCCC's.




Figure 8: 40X of a microfluidic channel showing captured NSCLC MCCCs from a Stage IV NSCLC patient. Anti-EGFR-AF594 (red), anti-CD44-FITC (green), Hoechst 33342 stained nuclei (blue), SMART-Coating™ (purple).




Figure 9: Gallery of MCCCs captured from patient whole blood samples. Heterogeneous shapes, sizes, and cell composition of the clusters are observed. The design of the chip and properties of the SMART-Coating™ allows culturing of viable MCCCs. Furthermore, we demonstrate efficient immunostaining procedures that enable advanced confocal imaging for detailed phenotypic characterization, and *ex vivo* drug testing.

RNA Seq on patient-captured MCCCs

RNA was extracted from sample INMC-44.6 (*i.e.* sample LC6). This sample corresponded to a 64-years old female patient, diagnosed by transthoracic biopsy with a TTF-1 positive stage IV lung adenocarcinoma that had metastasized to the brain. Three RNA aliquots of the LC6 sample were processed for library preparation using SMART-Seq® HT Kit (Takara) and Illumina Nextera reagents for fragmentation. RNA sequencing was performed at The Scripps Research Institute NGS core facility. Resulting FASTQ files were downloaded and imported into R. All subsequent analyses were performed using R/Bioconductor packages. Other datasets used in addition were: GSE51827: a study comparing gene expression profiles of clusters (CCLs) and single cells (SCs) in breast cancer.^{5,11} GSE74639: RNA-Seq of lung tumor circulating single cells (SCs) and primary lung tumor cells (PTs) from an orthotopic lung xenograft model.¹²




Figure 10: A) The overall count distribution was similar for the three replicates LC6.4, LC6.3, and LC6.6. B) INMC-44.6 is positive for the adenocarcinoma marker TTF-1. We confirmed expression of the gene "NKX2.1" which encodes the TTF-1 protein in MCCCs (together with CTCs and PTs). C) Heatmap showing the relationship of MCCCs to lung CTCs and primary lung tumor cells (PT), using genes that are differentially expressed between CTCs and PTs.

Detailed analysis of the RNA seq data obtained from capture MCCCs confirmed tissue of origin (not shown) as well as unique expression patterns specific to the patient (TTF-1). Furthermore, differentially expressed gene profiles unique to the MCCCs provides evidence that these vectors of metastasis require closer observation. Recent work by Aceto, *et al.*¹³, has highlighted highly promising therapeutic value in targeting MCCCs directly in the fight against metastasis.

Acknowledgements

We are truly grateful for the participation of the patients in our study and the support of Drs. Sandip Patel and Greg Botta from UCSD Moores Cancer Center. Thank you also to Steven Head, Director of the NGS Core at TSRI, for aid with RNA seq sample preparation, and James Evans from PhenoVista for support and encouragements. This work was supported by the NIH/NCI grant 1R43CA261362.

References

- Wang, *et al.*, Survival Mechanisms and Influence Factors of Circulating Tumor Cells. *BioMed Res. Int.*, 2018
- Rodrigues *et al.*, Circulating Tumor Cells: Come Together, Right Now, Over Metastasis. *Cancer Discovery*, 2019
- Rigg, *et al.*, Cancer Metastasis: A Reappraisal of Its Underlying Mechanisms and Their Relevance to Treatment. *Annual Review of Pathology*, 2018
- Lim, *et al.*, Circulating Tumor Cell Clusters Are Cloaked with Platelets and Correlate with Poor Prognosis in Unresectable Pancreatic Cancer. *Cancers*, 2021
- Aceto, *et al.*, Circulating Tumor Cell Clusters Are Oligoclonal Precursors of Breast Cancer Metastasis. *Cell*, 2014
- Costa, *et al.*, Analysis of a Real-World Cohort of Metastatic Breast Cancer Patients Shows Circulating Tumor Cell Clusters (CTC-clusters) as Predictors of Patient Outcomes. *Cancers*, 2020
- Murthi, *et al.*, Poor Prognosis Indicated by Venous Circulating Tumor Cell Clusters in Early-Stage Lung Cancer. *Cancer Res.*, 2017
- Sharma, *et al.*, Heterotypic Clustering of Circulating Tumor Cells and Circulating Cancer-associated Fibroblasts Facilitates Breast Cancer Metastasis. *Breast Can. Res. Treatment*, 2021
- Shaner *et al.*, Design and production of a novel microfluidic device for the capture and isolation of circulating tumor cell clusters. *AIIP Advances*, 2019
- UCSD Moores Cancer Center, IRB 181755, Project ID# INMC-044
- Yu, *et al.*, *Ex vivo* culture of circulating breast tumor cells for individualized testing of drug susceptibility. *Can. Therapy.*, 2014
- Zheng, *et al.*, Expression of β-globin by cancer cells promotes cell survival during blood-borne dissemination. *Nat. Com.*, 2017
- Gkoutis, *et al.*, Circulating Tumor Cell Clustering Shapes DNA Methylation. *Cell*, 2019